



Description: *GMO Chek RUR Soya Kit*

Part Number: 7100000

Lot Number: 10A1004

Expiration Date: February 2011

Storage Condition: Refrigerated, 2 to 8 °C

Quality Control Testing

Test	Results	Pass/Fail
Mean OD of 0% Reference Standard	0.079	Pass
Mean OD of 2.5% Reference Standard	1.2	Pass
Correlation Coefficient of Reference Standards	0.999	Pass
Mean OD of 0 ng/mL Protein Standard	0.112	Pass

The values above are the mean of all Quality Control results obtained for this lot. Results are affected by temperature, time, and pipette accuracy. Test performed in a controlled temperature setting with calibrated instrumentation. Individual results may vary from results listed above.

Constance Viering

Director of Kit Operations

Insert
Approval
Stamp Here

Alan McQuillin

Director of Quality

