



Certificate of Analysis

Description: Trait Chek RUR Lateral Flow Test Strip
Part Number: 3000082
Lot Number: 09J1009X
Expiration Date: March 2011

Storage Condition: Room Temperature

Quality Control Testing

Test	Results	Specification	Pass/Fail
Control line visible on all strips	Visible on all strips	Visible on all strips	Pass
Negative Control @ 5 Minutes (Non-GMO Soybean)	No Test Line present	No Test Line present	Pass
Positive Control @ 5 Minutes (RUR Soybean)	Test Line visible	Test Line visible	Pass
Negative Control @ 60 Minutes (Non-GMO Soybean)	No Test Line present	No Test Line present	Pass

The values above are the mean of all Quality Control results obtained for this lot.

Constance Viering
Director of Kit Operations

Insert
Approval
Stamp Here

9/21/2009

Alan McQuillin
Director of Quality