



**Description:** *Rapid Chek E. coli Media*

**Part Number:** 7000161, 7000165

**Lot Number:** C10-50-2

**Expiration Date:** April 2014

**Storage Condition:** Room Temperature

**Physical and Chemical Characteristics**

**Appearance (Dry):** Homogeneous, free flowing w/lumps, yellow beige

**Color/Clarity (Reconstituted):** Yellow/clear

**pH (2.5% solution at 25°C: 7.0 +/- 0.2):** 7.1

**Quality Control Testing**

Test	Results	Specification	Pass/Fail
Escherichia coli ATCC 45750 Control Organism	Good Growth	Good Growth	Pass
Escherichia coli ATCC 35150 Control Organism	Good Growth	Good Growth	Pass
Bacillus cereus ATCC 11778 Control Organism	Inhibited Growth	Inhibited Growth	Pass
Proteus vulgaris ATCC 8427 Control Organism	Inhibited Growth	Inhibited Growth	Pass
Helococcus kunzii ATTC 152 Control Organism	Inhibited Growth	Inhibited Growth	Pass
Streptococcus bovis ATCC 9809 Control Organism	Partial Inhibition	Partial Inhibition	Pass

Constance Viering

Director of Kit Operations

Insert Approval Stamp Here

Alan McQuillin

Director of Quality