



## Certificate of Analysis

**Description: Rapid Chek E. coli Media**

**Part Number: 7000161, 7000165**

**Lot Number: G09-13**

**Expiration Date: August 2013**

**Storage Condition: Room Temperature**

### Physical and Chemical Characteristics

**Appearance (Dry):** Homogeneous, free flowing w/lumps, yellow beige

**Color/Clarity (Reconstituted):** Yellow/clear

**pH (2.5% solution at 25°C: 7.0 +/- 0.2):** 6.96

### Quality Control Testing

Test	Results	Specification	Pass/Fail
<b>Escherichia coli ATCC 45750 Control Organism</b>	Good Growth	Good Growth	Pass
<b>Escherichia coli ATCC 35150 Control Organism</b>	Good Growth	Good Growth	Pass
<b>Bacillus cereus ATCC 11778 Control Organism</b>	Inhibited Growth	Inhibited Growth	Pass
<b>Proteus vulgaris ATCC 8427 Control Organism</b>	Inhibited Growth	Inhibited Growth	Pass
<b>Helococcus kunzii ATCC 152 Control Organism</b>	Inhibited Growth	Inhibited Growth	Pass
<b>Streptococcus bovis ATCC 9809 Control Organism</b>	Partial Inhibition	Partial Inhibition	Pass

Constance Viering  
Director of Kit Operations

Insert  
Approval  
Stamp Here

8/14/2009

Alan McQuillin  
Director of Quality