



Certificate of Analysis

Description: Rapid Chek E. coli Media

Part Number: 7000161, 7000165

Lot Number: H09-41-2

Expiration Date: September 2013

Storage Condition: Room Temperature

Physical and Chemical Characteristics

Appearance (Dry): Homogeneous, free flowing w/lumps, yellow beige

Color/Clarity (Reconstituted): Yellow/clear

pH (2.5% solution at 25°C: 7.0 +/- 0.2): 6.95

Quality Control Testing

Test	Results	Specification	Pass/Fail
Escherichia coli ATCC 45750 Control Organism	Good Growth	Good Growth	Pass
Escherichia coli ATCC 35150 Control Organism	Good Growth	Good Growth	Pass
Bacillus cereus ATCC 11778 Control Organism	Inhibited Growth	Inhibited Growth	Pass
Proteus vulgaris ATCC 8427 Control Organism	Inhibited Growth	Inhibited Growth	Pass
Helococcus kunzii ATCC 152 Control Organism	Inhibited Growth	Inhibited Growth	Pass
Streptococcus bovis ATCC 9809 Control Organism	Partial Inhibition	Partial Inhibition	Pass

Constance Viering
Director of Kit Operations

Insert
Approval
Stamp Here

8/14/2009

Alan McQuillin
Director of Quality