



## Certificate of Analysis

**Description: Rapid Chek E. coli Media**

**Part Number: 7000161, 7000165**

**Lot Number: H08-41**

**Expiration Date: September 2012**

**Storage Condition: Room Temperature**

### Physical and Chemical Characteristics

**Appearance (Dry):** Homogeneous, free flowing w/lumps, yellow beige

**Color/Clarity (Reconstituted):** Yellow/clear

**pH (2.5% solution at 25°C: 7.0 +/- 0.2):** 6.95

### Quality Control Testing

Test	Results	Specification	Pass/Fail
<b>E.coli ATCC 45750 Control Organism</b>	Good Growth	Good Growth	Pass
<b>E.coli ATCC 35150 Control Organism</b>	Good Growth	Good Growth	Pass
<b>B.cereus ATCC 11778 Control Organism</b>	Inhibited Growth	Inhibited Growth	Pass
<b>P.vulgaris ATCC 8427 Control Organism</b>	Inhibited Growth	Inhibited Growth	Pass

Constance Viering  
Director of Kit Operations

Insert  
Approval  
Stamp Here

3/25/2009

Alan McQuillin  
Director of Quality