



Description: *RapidChek® Base Enrichment Media for Listeria*

Part Number: 7000172, 7000177

Lot Number: C10-59

Expiration Date: April 2014

Storage Condition: Room Temperature

Physical and Chemical Characteristics

Appearance (Dry): Homogeneous, free flowing, yellow beige to tan powder

Color/Clarity (Reconstituted): Yellow to amber/clear to slight haze

pH (5.3% solution at 25°C: 7.2 +/- 0.2): 7.2

Quality Control Testing

Test	Results	Specification	Pass/Fail
L. innocua ATCC 33090	Good Growth	Good Growth	Pass
L. monocytogenes ATCC 7644	Good Growth	Good Growth	Pass
L. monocytogenes ATCC 15313	Good Growth	Good Growth	Pass
L. monocytogenes ATCC 19112	Good Growth	Good Growth	Pass
L. monocytogenes ATCC 19115	Good Growth	Good Growth	Pass
E. faecalis ATCC 10741	No Growth	No Growth	Pass
E. coli ATCC 25922	No Growth	No Growth	Pass
P.aeruginosa ATCC 27853	No Growth	No Growth	Pass
S.aureus ATCC 25923	No Growth	No Growth	Pass

Constance Viering
Director of Kit Operations



Alan McQuillin
Director of Quality