



Certificate of Analysis

Description: Rapid Chek Listeria Media 5.3kg
Part Number: 7000177
Lot Number: K07-31
Expiration Date: December 2011
Storage Condition: Room Temperature

Physical and Chemical Characteristics

Appearance (Dry): Homogeneous, free flowing, beige
Color/Clarity (Reconstituted): yellow/clear
pH (5.3% solution at 25°C): 7.16

Quality Control Testing

Test	Results	Specification	Pass/Fail
L.monocytogenes ATCC 7644 Control Organism	Good Growth	Good Growth	Pass
L.monocytogenes ATCC 19115 Control Organism	Good Growth	Good Growth	Pass
L.monocytogenes ATCC 19112 Control Organism	Good Growth	Good Growth	Pass
S.aureus ATCC 25923	No Growth	No Growth	Pass
P.aeruginosa ATCC 27853	No Growth	No Growth	Pass
E. coli ATCC 25922	No Growth	No Growth	Pass

This product has been manufactured, processed and packaged in accordance with Quality Systems Regulation, 21 CFR Part 820. This lot meets the criteria for approval as defined for this product. It should be stable to the expiration date on the label when stored according to package specifications.

Quality Assurance

111 Pencader Drive | Newark, DE 19702 | USA
(800) 481-9737 | +1(302) 607-3080