



Certificate of Analysis

Product Description: RapidChek® Listeria Test Strip

Part Numbers: 3000019

Lot Number: 09B1029X

Expiration Date: August 2010


Storage Condition: Room Temperature

Quality Control Testing

Test	Results	Specification	Pass/Fail
Control line visible on all strips	Visible on all strips	Visible on all strips	Pass
Negative Control @ 10 Minutes	No Test Line present	No Test Line present	Pass
Negative Control @ 60 Minutes	No Test Line present	No Test Line present	Pass
Positive Control @ 10 Minutes	Test Line visible	Test Line visible	Pass

The values above are the mean of all Quality Control results obtained for this lot

Constance Viering
Director of Kit Operations



3/3/2009

Alan McQuillin
Director of Quality