



Certificate of Analysis

Description: Rapid Chek Medium for Salmonella

Part Number: 7000167

Lot Number: E08-02-01

Expiration Date: June 2012

Storage Condition: Room Temperature

Physical Characteristics

Appearance (Dry): Homogeneous, free flowing w/lumps, light beige

Color/Clarity: Light yellow/clear

pH (2.0% Solution at 25°C: 7.0 +/- 0.1): 7.00

Quality Control Testing

Test	Results	Specification	Pass/Fail
S.enteritidis ATCC 13076 Control Organism	Good Growth	Good Growth	Pass
S.montevidaea SDI ARS32 Control Organism	Good Growth	Good Growth	Pass
S.typhimurium ATCC 14028 Control Organism	Good Growth	Good Growth	Pass

Quality Assurance

111 Pencader Drive | Newark, DE 19702 | USA

(800) 481-9737 | +1(302) 607-3080