



## Certificate of Analysis

**Description: Rapid Chek Media for Salmonella**

**Part Number: 7000167**

**Lot Number: G09-34-3**

**Expiration Date: August 2013**

**Storage Condition: Room Temperature**

**Physical Characteristics**

**Appearance (Dry):** Homogeneous, free flowing w/lumps, light beige

**Color/Clarity:** Light yellow/clear

**pH (2.0% Solution at 25°C: 7.0 +/- 0.1):** 7.0

**Quality Control Testing**

Test	Results	Specification	Pass/Fail
<b>S.enteritidis ATCC 13076 Control Organism</b>	Good Growth	Good Growth	Pass
<b>S.montevidea SDI ARS32 Control Organism</b>	Good Growth	Good Growth	Pass
<b>S.typhimurium ATCC 14028 Control Organism</b>	Good Growth	Good Growth	Pass

Constance Viering  
Director of Kit Operations

Insert  
Approval  
Stamp Here

2/3/2010

Alan McQuillin  
Director of Quality