



# Certificate of Analysis

**Description:** Rapid Chek Medium for Salmonella

**Part Number:** 7000167

**Lot Number:** L07-15

**Expiration Date:** January 2012

**Storage Condition:** Room Temperature

## Physical Characteristics

**Appearance (Dry):** Homogeneous, free flowing w/lumps, light beige

**Color/Clarity:** Light yellow/clear

**pH (2.0% Solution at 25°C):** 6.92

## Quality Control Testing

Test	Results	Specification	Pass/Fail
S.enteritidis ATCC 13076 Control Organism	Good Growth	Good Growth	Pass
S.montevidaea SDI ARS32 Control Organism	Good Growth	Good Growth	Pass
S.typhimurium ATCC 14028 Control Organism	Good Growth	Good Growth	Pass

## Quality Assurance

111 Pencader Drive | Newark, DE 19702 | USA

(800) 481-9737 | +1(302) 607-3080