



**Description:** *Rapid Chek Salmonella Primary Enrichment Medium*

**Part Number:** 7000184, 7000192

**Lot Number:** C10-54-1

**Expiration Date:** April 2014

**Storage Condition:** Room Temperature

***Physical and Chemical Characteristics***

Appearance (Dry): Homogeneous, free flowing with soft lumps, light beige

Color/Clarity (Reconstituted): Light yellow/clear

pH (2.0% solution at 25°C: 7.0 +/- 0.1): 7.1

***Quality Control Testing***

Test	Results	Specification	Pass/Fail
Salmonella enteritidis ATCC 13076 Control Organism	Good Growth	Good Growth	Pass
Salmonella montevideo ARS32 Control Organism	Good Growth	Good Growth	Pass
Salmonella typhimurium ATCC 14028 Control Organism	Good Growth	Good Growth	Pass

Constance Viering

Director of Kit Operations

Insert  
Approval  
Stamp Here

Alan McQuillin

Director of Quality