



Description: *Rapid Chek Salmonella Primary Enrichment Medium*

Part Number: 7000184, 7000192

Lot Number: D09-05-2

Expiration Date: May 2013

Storage Condition: Room Temperature

Physical and Chemical Characteristics

Appearance (Dry): Homogeneous, free flowing with soft lumps, light beige

Color/Clarity (Reconstituted): Light yellow/clear

pH (2.0% solution at 25°C: 7.0 +/- 0.1): 7.0

Quality Control Testing

Test	Results	Specification	Pass/Fail
Salmonella enteritidis ATCC 13076 Control Organism	Good Growth	Good Growth	Pass
Salmonella montevideo ARS32 Control Organism	Good Growth	Good Growth	Pass
Salmonella typhimurium ATCC 14028 Control Organism	Good Growth	Good Growth	Pass

Constance Viering

Director of Kit Operations

Insert
Approval
Stamp Here

Alan McQuillin

Director of Quality