



Certificate of Analysis

Description: *RapidChek® SELECT™ Salmonella Primary Enrichment Media*

Part Number: 7000184, 7000192

Lot Number: F10-18

Expiration Date: July 2014

Storage Condition: Room Temperature

Physical and Chemical Characteristics

Appearance (Dry): Homogeneous, free flowing w/ soft lumps, light beige

Color/Clarity (Reconstituted): Light yellow to tan or amber/clear

pH (2.0% solution at 25°C: 7.0 +/- 0.2): 7.0

Quality Control Testing

Test	Results	Specification	Pass/Fail
Salmonella enteritidis ATCC 13076	Good Growth	Good Growth	Pass
Salmonella montivideo SDI ARS32	Good Growth	Good Growth	Pass
Salmonella typhimurium ATCC 14028	Good Growth	Good Growth	Pass

Constance Viering

Director of Kit Operations



Alan McQuillin

Director of Quality

111 Pencader Drive | Newark, DE 19702 USA | +1 (800) 481-9737 | +1 (302) 607-3080

www.sdix.com