



## Certificate of Analysis

---

**Description: Microtox 2% Diluent - 1 Liter**

**Part Number: AZF686011**

**Lot Number: 09M1022**

**Expiration Date: December 2012**

**Storage Condition: Room Temperature**

### Quality Control Testing

Test	Results	Specification	Pass/Fail
Comparison % - 5 Min	1.0%	$\leq \pm 5\%$	Pass
Comparison % - 15 Min	1.1%	$\leq \pm 5\%$	Pass
Conductivity result	2.16%	1.68% - 2.52%	Pass

The values above are the mean of all Quality Control results obtained for this lot.

All comparison results are obtained using the Quality Control Procedure with an approved lot of Microtox 2% Diluent on the Model 500.

Constance Viering  
Director of Kit Operations

Insert  
Approval  
Stamp Here

9/8/2010

Alan McQuillin  
Director of Quality