

Description: Microtox 0.01% Reconstitution Solution**Part Number:** AZF686016**Lot Number:** 10H1043**Expiration Date:** August 2013**Storage Condition:** Room Temperature**Quality Control Testing**

Test	Results	Specification	Pass/Fail
Comparison % - 5 Min	0.3%	$\leq \pm 5\%$	Pass
Comparison % - 15 Min	1.2%	$\leq \pm 5\%$	Pass
Conductivity result	0.010%	0.008% - 0.014%	Pass

The values above are the mean of all Quality Control results obtained for this lot.

All comparison results are obtained using the Quality Control Procedure with an approved lot of Microtox 0.01% Diluent on the Model 500.



Constance Vierung

Director of Kit Operations



Alan McQuillin

Director of Quality