

Anti- human PAX5 Mouse Monoclonal Antibody

Clone: 3A7 REF RU00084

Intended use

Anti- human PAX5 (Clone: 3A7) Mouse Monoclonal Primary Antibody is intended for research use only. Not for use in diagnostic procedures. Not for human or animal consumption.

Background

Paired box 5 is a member of the paired box (PAX) family of transcription factors. The central feature of this gene family is a novel, highly conserved DNA-binding motif, known as the paired box. PAX proteins are important regulators in early development, and alterations in the expression of their genes are thought to contribute to neoplastic transformation. This gene encodes the B-cell lineage specific activator protein that is expressed at early, but not late stages of B-cell differentiation. Its expression has also been detected in developing CNS and testis and so the encoded protein may also play a role in neural development and spermatogenesis. This gene is located at 9p13, which is involved in t(9;14)(p13;q32) translocations recurring in small lymphocytic lymphomas of the plasmacytoid subtype, and in derived large-cell lymphomas. This translocation brings the potent E- μ enhancer of the IgH gene into close proximity of the PAX5 promoter, suggesting that the deregulation of transcription of this gene contributes to the pathogenesis of these lymphomas. Alternatively spliced transcript variants encoding different isoforms have been described but their biological validity has not been determined. [provided by RefSeq, Jul 2008].

Alternative names: ALL3; BSAP

Reagent provided

Anti-human PAX5 Mouse Monoclonal Primary Antibody (Clone: 3A7) is provided in liquid form in 20mM Sodium phosphate, 150mM Sodium chloride, 0.2% BSA, 0.09% Sodium azide, pH 7.4. The isotype of the antibody is IgG1,k. The protein concentration is approximately 0.6 +/- 0.05 mg/mL. For immunohistochemistry, the primary antibody may be used at a working dilution of 1:100 – 1:200 for formalin-fixed, paraffin-embedded human tissue. It can be dependent upon the detection system used. These are guidelines only, and optimal dilutions should be determined by the individual laboratory.

Immunogen

Full length human recombinant protein of human PAX5 (NP_057953) was produced in E. coli.

Specificity

The specificity of the anti-human PAX5 Mouse Monoclonal Primary Antibody was established on known human tonsil tissue. The anti-human PAX5 presented no staining on human brain tissue and positive staining on human tonsil tissue using immunohistochemical (IHC) test methods.

Precautions

1. This product contains sodium azide (NaN_3), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN_3 may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
2. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
3. Unused reagents should be disposed of according to local, State, and Federal regulations.
4. Suitability for specific application may vary and it is the responsibility of the end user to determine the appropriate application for use and stability.

Storage

Store at 2-8°C. If reagents are stored under conditions other than those specified in the package insert, they must be verified by the user.

Your Antibody Partner – Innovating to a Healthier Tomorrow

SDIX, LLC • 111 Pencader Drive • Newark, DE 19702 • 800-544-8881 • 302-456-6789 • www.sdix.com • Rev.01/22